



Catalog

Human Biological Products and Services

WHOLE BLOOD



SERUM & PLASMA



DISEASE-STATE



WWW.EQUITECH-ENTERPRISES.COM | HUMANSALES@EQUITECH-BIO.COM

Contents

About Us	2
Ordering Information	2
Human IgG and Albumin	3
Human Immunoglobulin G (Human IgG: Lyophilized Powder)	3
Human Immunoglobulin G (Human IgG: Liquid, Sterile-Filtered)	3
Human Serum Albumin: Lyophilized Powder	3
Human Whole Blood	3
Human Platelet Concentrate	3
Whole Blood from Designated Individuals or Disease-State Donors	3
Human Umbilical Cord Whole Blood Pooled & Individual Units	3
Human Serum Products	4
Serum Off The Clot (Normal Donors) Pooled & Individual	4
Serum AB Off The Clot Pooled & Individual	4
Serum AB Off The Clot—Male only Pooled	4
Serum Off The Clot with Human Anti-Mouse Antibodies (HAMA)	5
Dextran-Coated Charcoal-Stripped Human Serum Pooled	5
Serum from Designated Individuals or Disease-State Donors Pooled	5
Umbilical Cord Serum Pooled & Individual Units	5
Human Plasma Products	6
Human Plasma Defibrinated (Df), Delipidized (Dp) Pooled	6
Human Plasma Pooled & Individual	6
Plasma from Designated Individuals or Disease-State Donors (any anti-coagulant)	7
Unfiltered Plasma from Designated Individuals or Disease-State Donors	7

ABOUT US

Our Blood Collection Center

The EEI Blood Collection Center is FDA Registered, and provides Informed Consent Collections. Our Staff Donor Recruiters ensure a diversified donor population that offers a wide range of age, ethnicity, and disease state. All biologicals collected receive viral screening including HIV, Hepatitis B & C, Syphilis, and other selected markers.

Data and Demographics

All volunteers receive, and are qualified by, a fully comprehensive Donor Health Questionnaire from which data and demographics information is offered upon request. EEI can provide custom specimen collections or bulk quantities ranging from a single unit to multiple liters.

Our Quality Policy

It is Equitech Enterprise's goal and mission to be recognized as a leading supplier of quality blood, and blood components to the research industry. To achieve this goal, it is our policy to set and achieve quality standards that allow us to provide exceptional material that is specific to your needs. Our services begin the moment you call us or email us at humansales@equitech-bio.com where you will receive answers to your questions, discuss your requirements. We will provide assurance that you'll receive the product you need, when you need it, and at a reasonable cost.

Aliquots, Pools, and Containers

EEI will collect materials in aliquot or pool into containers to your specifications. Plastic bottle, plastic bag, vials, and tubes.

Storage, Shipping, and Delivery

Materials can be frozen at minus 20 degrees C or kept refrigerated between 2 and 8 degrees C, or stored at room temperature where appropriate. Products will be shipped on dry ice or cold packs, or sent ambient temperature as requested. Overnight freight charges can be charged directly to your Federal Express, or courier of choice account.

ORDERING INFORMATION

Please direct orders to:

Equitech Enterprises, Inc.
1436-A Sidney Baker Street,
Kerrville, Texas 78028
E-mail: humansales@equitech-bio.com
Website: www.equitech-enterprises.com

Please direct inquiries to:

Email: humansales@equitech-bio.com

When placing an order, please specify customer account number, shipping address, billing address, catalog number, quantity and size, and your purchase order number. Please send a confirming purchase order for all email requests, and indicate clearly that it is for confirmation. For information on blanket orders, standing orders, prepaid accounts and other special ordering arrangements, or larger quantities, please call us. Catalog prices are for small quantities only.

TERMS OF PAYMENT AND DELIVERY

ORDERS

All orders and agreements are binding only if confirmed in writing at our Kerrville, Texas office or through delivery of merchandise. When ordering, please be certain to indicate both catalog number and description of product. **Equitech Enterprises, Inc. has a minimum of \$100 for domestic orders, and \$250 for international orders.** Orders below this value will have an administrative charge added.

SHIPMENTS

We recommend that customers give specific instructions regarding shipment of their order. In the absence of instructions we will ship the most practical way under conditions existing at the time of shipment.

TERMS & PRICES

Payment is due within 30 days of date of invoice with the exception that for international orders special payment terms need to be established. Prices are subject to change without notice. All orders are subject to final acceptance by Equitech Enterprises, Inc. at the prices in effect at the time of shipment. Applicable estimate or quote number must be referenced at the time an order is placed in order to receive quoted special pricing. Prices are net FOB Kerrville, Texas plus shipping and handling.

CLAIMS FOR DAMAGE IN TRANSIT

Our shipping responsibility ceases with safe delivery in good condition to the transportation company. Claims by the customer for damage in transit should be made to the transportation company as follows: If upon delivery the outside of the packaging case shows evidence of rough handling or damage, the transportation company's agent should be requested to make a "Received in Bad Order" notation on the delivery receipt. If there is no exterior evidence of rough handling upon delivery but concealed damage is evident upon unpacking the shipment, the transportation company should, within 48 hours of delivery, be requested to make out a "Bad Order" report. This procedure is necessary in order to maintain the right of recovery from the carrier. In addition, Equitech Enterprises, Inc. should be informed within 10 days of receipt of damaged goods. Non-receipt claims on lost shipment should be made within 30 days of original order. Insurance is available for shipments up to \$6,000.

CLAIMS FOR SHORTAGES

We use extreme care in selection, checking, rechecking and packing to eliminate the possibility of error. If any shipment discrepancies are discovered, claims must be made to the place of shipment within 10 days to be considered. Please check all packing material thoroughly upon receipt of shipment.

RETURNING MATERIAL

No product may be returned without authorization from Equitech Enterprises, Inc. This authorization is necessary to ensure correct return of material and issuance of credit or exchange of material. Requests for permission to return products must be made within ten days of receipt. Returned products should be clearly marked with sender's name and address, carefully packed, and shipped prepaid. Products must be returned, unopened, in the same condition as received and within 15 days of original shipment by Equitech Enterprises, Inc.

HUMAN IgG AND ALBUMIN

The human IgG and albumin shown below are only available from our sister company Equitech-Bio, Inc. You can contact them on 1-800-259-0591 or email sales@equitech-bio.com, or from their website: www.equitech-bio.com.

Human Immunoglobulin G (Human IgG): Lyophilized Powder

Catalog No.	Size
SLH56-0001	1 g
SLH56-0010	10 g
SLH56-0100	100 g

Human Immunoglobulin G (Human IgG): Liquid, Sterile Filtered

Catalog No.	Size
H60-0001	1 g
H60-0010	10 g
H60-0100	100 g

Human Serum Albumin: Lyophilized Powder

- Prepared from human serum
- Purity: ≥ 96%
- Protein level: ≥ 96%
- IgG level: not detected

Catalog No.	Size
HSA62-0010	10 g
HSA62-0100	100 g
HSA62-1000	1000 g

WHOLE BLOOD PRODUCTS

Human Platelet Concentrate

Catalog No.	Size
HPT-0001	1ml

Whole Blood from Designated Individuals or Disease-State Donors

Based on anticoagulants and volume

- "As is" from donor
- Storage conditions: -4°C ± 2°C
- Each unit is tested for Anti-HIV 1/2, Anti-HCV, HBs Ag and RPR.

Catalog No.	Size
HWB-Unit	Donor Unit

Human Umbilical Cord Whole Blood

Based on anticoagulants and volume

Mothers tested negative for HIV 1/2, HBs Ag, Anti-HCV and RPR. Collected in the USA from hospitals with signed donor approval. Various levels of alpha-fetoprotein (AFP) are available; please inquire.

Pooled

Catalog No.	Size
HUCWB	Varies

Individual Units

Catalog No.	Size
IHUCWB	Varies

HUMAN SERUM PRODUCTS

Equitech Enterprise's human serum and plasma products are collected in donor centers located in the United States. Each unit is tested for Anti-HIV 1/2, Anti-HCV, HBs Ag and RPR.

Human Serum Off The Clot (Normal Donors)

Sterile-Filtered Serum

- Filtered to 0.2 micron
- Sterility tested x 14 days
- Low hemoglobin levels
- Storage conditions: -20°C ± 5°C

Pooled

Catalog Number	Size
SHS-0050	50ml
SHS-0100	100 ml
SHS-0500	500 ml
SHS-1000	1 Liter

Individual Units

Catalog Number	Size
ISHS-0050	50ml
ISHS-0100	100 ml
ISHS-0500	500 ml
ISHS-1000	1 Liter

Unfiltered Serum

- "As is" from patient
- Low hemoglobin levels
- Storage conditions: -20°C ± 5°C
- Tested for Anti-HIV 1/2, Anti-HCV, HBs Ag and RPR.

Pooled

Catalog Number	Size
HS-0050	50 ml
HS-0100	100 ml
HS-0500	500 ml
HS-1000	1 Liter

Individual Units

Catalog Number	Size
IHS-0050	50 ml
IHS-0100	100 ml
IHS-0500	500 ml
IHS-1000	1 Liter

Human Serum AB Off The Clot

Sterile-Filtered Serum

- Filtered to 0.2 micron
- Sterility tested x 14 days
- Low hemoglobin levels
- Storage conditions: -20°C ± 5°C

Pooled

Catalog No.	Size
SHSAB-0100	100 ml
SHSAB-0500	500 ml

Individual Units

Catalog No.	Size
ISHSAB-0100	100 ml
ISHSAB-0500	500 ml

Unfiltered Serum

- "As is" from donor
- Low hemoglobin levels
- Storage conditions: -20°C ± 5°C
- Tested for Anti-HIV 1/2, Anti-HCV, HBs Ag and RPR.

Pooled

Catalog No.	Size
HSAB-0100	100 ml
HSAB-0500	500 ml

Individual Units

Catalog No.	Size
IHSAB-0100	100 ml
IHSAB-0500	500 ml

Human Serum AB Off The Clot - Male Only

Sterile-Filtered Serum

- Filtered to 0.2 micron
- Sterility tested x 14 days
- Low hemoglobin levels
- Storage conditions: -20°C ± 5°C

Pooled

Catalog No.	Size
SHSABM-0100	100 ml
SHSABM-0500	500 ml

Unfiltered Serum

- "As is" from donor
- Low hemoglobin levels
- Storage conditions: -20°C ± 5°C
- Tested for Anti-HIV 1/2, Anti-HCV, HBs, Ag and RPR.

Pooled

Catalog No.	Size
HSABM-0100	100 ml
HSABM-0500	500 ml

Continued from page 4.

Individual Units

Catalog No.	Size
ISHSABM-0100	100 ml
ISHSABM-0500	500 ml

Individual Units

Catalog No.	Size
IHSABM-0100	100 ml
IHSABM-0500	500 ml

Human Serum Off The Clot with Human Anti-Mouse Antibodies (HAMA):

This serum has been tested for Human Anti-Mouse Antibodies. Different titers are available.

Individual Units Only

Catalog No.	Size
IHAMA-0100	1 ml

Dextran-Coated Charcoal-Stripped Human Serum

Charcoal Stripping reduces the concentration of steroid hormones in serum including as estradiol, progesterone, cortisol, testosterone T3 and T4.

- Sterile-filtered to 0.2 micron
- Sterility tested x 14 days
- Storage conditions: -20°C ± 5°C
- Each unit is tested for Anti-HIV 1/2, Anti-HCV, HBs Ag and RPR.

Pooled

Catalog No.	Size
SCHS-0100	100 ml
SCHS-0500	500 ml
SCHS-1000	1 Liter

Individual Units

Catalog No.	Size
ISCHS-0100	100 ml
ISCHS-0500	500 ml
ISCHS-1000	1 Liter

Serum from Designated Individuals or Disease-State Donors

Sterile-Filtered Serum

- Sterile-filtered to 0.2 micron
- Sterility tested x 14 days
- Low hemoglobin levels
- Storage conditions: -20°C ± 5°C

Unfiltered Serum

- "As is" from donor
- Low hemoglobin levels
- Storage conditions: -20°C ± 5°C
- Tested for Anti-HIV 1/2, Anti-HCV, HBs, Ag and RPR.

Sterile Filtered Pooled

Catalog No.	Size
SHS-0005 D	5 ml
SHS-0010 D	10 ml
SHS- Unit D	Unit

Unfiltered Pooled

Catalog No.	Size
HS-0005 D	5 ml
HS-0010 D	10 ml
HS-Unit D	Unit

Human Umbilical Cord Serum

Mothers tested negative for HIV 1/2, HBs Ag, Anti-HCV and RPR. Collected in the USA from hospitals with signed donor approval. Various levels of alpha-fetoprotein (AFP) are available; please inquire.

Equitech-Bio offers sera from donors with syphilis, Hepatitis B, Hepatitis C, arthritis, and elevated analytes for specific markers, and sensitivities for specific allergens.

Pooled

Catalog No.	Size
HUCS	Varies

Individual Units

Catalog No.	Size
IHUCS	Varies

HUMAN PLASMA PRODUCTS

Human Plasma Defibrinated (Df), Delipidized (Dp)

- Sterile-filtered to 0.2 micron
- Sterility tested x 14 days
- Storage conditions: -20°C ± 5°C
- Each unit is tested for Anti-HIV 1/2, Anti-HCV, HBs Ag and RPR.

Pooled

Catalog No.	Size
SDfDpHP-0100	100ml
SDfDpHP-0500	500ml

Individual Units

Catalog No.	Size
ISDfDpHP-0100	100ml
ISDfDpHP-0500	500ml

Charcoal-Stripped

Pooled

Catalog No.	Size
SCDfDpHP-0100	100ml
SCDfDpHP-0500	500ml

Charcoal-Stripped

Individual Units

Catalog No.	Size
ISCDfDpHP-0100	100ml
ISCDfDpHP-0500	500ml

Human Plasma

Sterile-Filtered Pooled

- Sterile-filtered to 0.2 micron
- Sterility tested x 14 days
- Storage conditions: -20°C ± 5°C

Catalog No.	Size
CITRATE	
SHPC-0100	100 ml
SHPC-0500	500 ml
SHPC-1000	1 Liter
CPD	
SHPCPD-0100	100 ml
SHPCPD-0500	500 ml
SHPCPD-1000	1 Liter
K2-EDTA	
SHPEK2-0100	100 ml
SHPEK2-0500	500 ml
SHPEK2-1000	1 Liter
K3-EDTA	
SHPEK3-0100	100 ml
SHPEK3-0500	500 ml
SHPEK3-1000	1 Liter
LITHIUM HEPARIN	LITHIUM HEPARIN
SHPLH-0100	100 ML
SHPLH-0500	500 ML
SHPLH-1000	1 Liter
SODIUM HEPARIN	SODIUM HEPARIN
SHPNH-0100	100 ML
SHPNH-0500	500 ML
SHPNH-1000	1 Liter

Unfiltered Pooled

- "As is" from donor
- Storage conditions: -20°C ± 5°C
- Tested for Anti-HIV 1/2, Anti-HCV, HBs Ag and RPR.

Catalog No.	Size
CITRATE	
HPC-0100	100 ml
HPC-0500	500 ml
HPC-1000	1 Liter
CPD	
HPCPD-0100	100 ml
HPCPD-0500	500 ml
HPCPD-1000	1 Liter
K2-EDTA	
HPEK2-0100	100 ml
HPEK2-0500	500 ml
HPEK2-1000	1 Liter
K3-EDTA	
HPEK3-0100	100 ml
HPEK3-0500	500 ml
HPEK3-1000	1 Liter
LITHIUM HEPARIN	LITHIUM HEPARIN
HPLH-0100	100 ML
HPLH-0500	500 ML
HPLH-1000	1 Liter
SODIUM HEPARIN	SODIUM HEPARIN
HPNH-0100	100 ML
HPNH-0500	500 ML
HPNH-1000	1 Liter

Sterile-Filtered Individual Units

Catalog No.	Size
CITRATE	
ISHPC-0100	100 ml
ISHPC-0500	500 ml
ISHPC-1000	1 Liter
CPD	
ISHPCPD-0100	100 ml
ISHPCPD-0500	500 ml
ISHPCPD-1000	1 Liter

Unfiltered Individual Units

Catalog No.	Size
CITRATE	
IHPC-0100	100 ml
IHPC-0500	500 ml
IHPC-1000	1 Liter
CPD	
IHPCPD-0100	100 ml
IHPCPD-0500	500 ml
IHPCPD-1000	1 Liter

Continued on page 7

Human Plasma Continued from p.6

Sterile-Filtered Individual Units

Catalog No.	Size
K2-EDTA	
ISHPEK2-0100	100 ml
ISHPEK2-0500	500 ml
ISHPEK2-1000	1 Liter
K3-EDTA	
ISHPEK3-0100	100 ml
ISHPEK3-0500	500 ml
ISHPEK3-1000	1 Liter
HEPARIN	
ISHPH-0100	100 ml
ISHPH-0500	500 ml
ISHPH-1000	1 Liter
LITHIUM HEPARIN	
ISHPLH-0100	100 ml
ISHPLH-0500	500 ml
ISHPLH-1000	1 Liter
SODIUM HEPARIN	
ISHPNH-0100	100 ml
ISHPNH-0500	500 ml
ISHPNH-1000	1 Liter

Unfiltered Individual Units

Catalog No.	Size
K2-EDTA	
IHPEK2-0100	100 ml
IHPEK2-0500	500 ml
IHPEK2-1000	1 Liter
K3-EDTA	
IHPEK3-0100	100 ml
IHPEK3-0500	500 ml
IHPEK3-1000	1 Liter
HEPARIN	
IHPH-0100	100 ml
IHPH-0500	500 ml
IHPH-1000	1 Liter
LITHIUM HEPARIN	
IHPLH-0100	100 ml
IHPLH-0500	500 ml
IHPLH-1000	1 Liter
SODIUM HEPARIN	
IHPNH-0100	100 ml
IHPNH-0500	500 ml
IHPNH-1000	1 Liter

Plasma from Designated Individuals or Disease State Donors—any anti-coagulant

Sterile-Filtered

- Sterility Tested x 14 days
- Sterile-filtered to 0.2 micron
- Storage conditions: -20°C ± 5°C
- Each unit is tested for Anti-HIV 1/2, Anti-HCV, HBs Ag and RPR.

Pooled

Catalog No.	Size
SHP - Unit D	Unit

Individual Units

Catalog No.	Size
ISHP - Unit D	Unit

Unfiltered Plasma from Designated Individuals or Disease-State Donors

Based on anticoagulants

- “As is” from donor
- Storage conditions: -20°C ± 5°C
- Each unit is tested for Anti-HIV 1/2, Anti-HCV, HBs Ag and RPR.

Pooled

Catalog No.	Size
HP-0005 D	5 ml
HP-0010 D	10 ml
HP-UNIT D	1 Unit

Individual Units

Catalog No.	Size
IHP-0005 D	5 ml
IHP-0010 D	10 ml
IHP-UNIT D	1 Unit



WWW.EQUITECH-ENTERPRISES.COM | HUMANSALES@EQUITECH-BIO.COM