



**BOVINE CALF SERUM
STERILE FILTERED
HEAT INACTIVATED**

TDS- SBCU32

Effective Date: 1/1/24

Rev. 000

Page 1 of 1

SPECIFICATIONS:

Source:	Bovine Calf Serum
Hemoglobin:	<40 mg/dl
Protein:	>6.0 g/dl
pH:	6.8 - 8.2
Osmolality:	260 – 360 mOsm/Kg
Filtration:	0.2 micron (Absolute)
Sterility:	No Growth

CATALOG NUMBER:

SBCU30-0100
SBCU32-0100
SBCU30-0500
SBCU32-0500

SIZES:

100ml
100ml Heat Inactivated
500ml
500ml Heat Inactivated

SHIPPING CONDITIONS: <0°C

STABILITY AND STORAGE:

This serum is stable for ≥ 5 years when stored $-20^{\circ} \text{C} \pm 5^{\circ}$. Avoid repeated freezing and thawing.

The raw materials for these products were collected from calves 12 months and younger in the USA in abattoirs registered with the USDA. The animals received ante and post-mortem inspections under a veterinarian's supervision and were apparently free from infectious and contagious diseases. At no time during collection or processing was the material commingled with any other material of animal origin.

This product is sold for use in vitro diagnostic, research, further manufacturing, or technical purposes only, and is not intended for direct use in diagnostic or therapeutic applications.