



**RABBIT SERUM
STERILE FILTERED**

TDS-SR30

Effective Date: March 8, 2023

Rev. 000

Page 1 of 1

SPECIFICATIONS:

Source:	Rabbit Serum, mixed sex, mixed breed
Hemoglobin:	<40 mg/dl
Protein:	≥5.0 g/dl
pH:	6.8 - 8.2
Osmolality:	260 – 360 mOsm/Kg
Filtration:	0.2 micron (Absolute)
Sterility:	No Growth

CATALOG NUMBER

SR30-0100

SR32-0100

SR30-0500

SR32-0500

SIZES

100ml

100ml Heat Inactivated

500ml

500ml Heat Inactivated

SHIPPING CONDITIONS: <0°C

STABILITY AND STORAGE:

This serum is stable for ≥ 3 years when stored -20° C ± 5°. Avoid repeated freezing and thawing.

The raw materials for these products were collected from healthy rabbits of French origin. The animals received ante and post-mortem inspections under a veterinarian's supervision and were apparently free from infectious and contagious diseases. At no time during collection or processing was the material commingled with any other material of animal origin. Product meets European Union requirements for production of technical blood products.

This product is sold for use in vitro diagnostic, research, further manufacturing, or technical purposes only, and is not intended for direct use in diagnostic or therapeutic applications.