

	<p align="center">GUINEA PIG SERUM STERILE FILTERED</p>	<p align="center">TDS-SGP30</p>
<p>Effective Date:</p>	<p>May 22, 2019</p>	<p>Rev. 000</p>
		<p>Page 1 of 1</p>

SPECIFICATIONS:

Source:	Guinea pigs, sexually mature, mixed sex, mixed breed
Hemoglobin:	<110 mg/dl
Protein:	>2.0 g/dl
pH:	6.8 - 8.2
Osmolality:	240 – 380 mOsm/Kg
Filtration:	0.2 micron (Absolute)
Sterility:	No Growth

CATALOG NUMBER

SGP30-0050

SGP32-0050

SGP30-0100

SGP32-0100

SGP30-0500

SGP32-0500

SIZES

50ml

50ml Heat Inactivated

100ml

100ml Heat Inactivated

500ml

500ml Heat Inactivated

SHIPPING CONDITIONS: <0°C

STABILITY AND STORAGE:

This serum is stable for ≥ 3 years when stored –20° C ± 5°. Avoid repeated freezing and thawing.

The material for this product was processed in the USA from the blood of healthy guinea pigs that had not been inoculated with or exposed to any infectious agents. At no time during collection or processing was the material commingled with any other material of animal origin.

This product is not to be used for direct or indirect exposure or inoculation into laboratory and domestic animals.

This product is sold for use in IVD, research, further manufacturing, or technical purposes only, and is not intended for direct use in diagnostic or therapeutic applications.